KAPS®

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Prothèse unicompartmentale de genou KAPS®

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EN

INSTRUCTIONS FOR USE

KAPS® unicompartmental knee prosthesis

Femoral component / Tibial insert / Tibial base / Tibial plateau / Cancellous screws

To be read carefully by the surgeon, theatre staff and all staff involved in the use and handling of the device.
The use of this device is strictly reserved to specialist practitioners with regular practice in knee replacement surgery and the implantation technique for unicompartmental femorotibial prostheses and must be carried out in premises that meet the asepsis conditions required and that comply with NF ISO 8828. Implants must be selected with particular attention to the user, taking account of the relevant indication, the patient's morphology and activity in order to minimise adverse effects and complications.
Any complications or the onset of side effects due to non-compliance with the recommendations mentioned in this guide and in the surgical technique, to material handling errors or an incorrect diagnosis would be entirely the responsibility of the surgeon.
X.NOV® does not give medical advice. Use of this product is strictly reserved to surgeons trained in knee surgery.

1 - Description

This box contains an implant (femoral component or tibial insert or tibial base or tibial plateau or cancellous screw) of the KAPS® unicompartmental knee prosthesis intended for the replacement of one of the tibiofemoral joints of the knee in a unicompartmental knee prosthesis operation.
X.NOV® knee prosthesis implants must be used only in combination with the other KAPS® knee prosthesis elements as intended by the manufacturer in this leaflet or in the surgical technique.
The components can be fixed with or without bone cement according to the information on the outer label of the box.
This device is intended for single use only.
Re-implantation of an explanted device is forbidden.

2 - Material

The metallic components are made from a cobalt chrome molybdenum alloy in compliance with ISO 5832-4.
The tibial insert as well as the tibial base to be cemented are made of UHMWPE polyethylene, in compliance with ISO 5834-1 and ISO 5834-2.
The osteoconductive coating of implants to be fitted without bone cement is made of a layer of pure titanium in compliance with ASTM F1580 and ISO 13179-1 and a layer of calcium hydroxyapatite in compliance with ASTM F1185 and ISO 13779-1 and 13779-2.
The cancellous bone screws are made of stainless steel in compliance with ISO 5832-9.
Bone cement: refer to the manufacturer's leaflet.

3 - Intended use

Unicompartmental arthroplasty of the knee is intended to improve the mobility of the patient and reduce pain by replacing one of the damaged femorotibial joints in patients with sufficient bone stock to attach and support the components of the prosthesis.

4 - Indications

- Destruction of the internal or external tibiofemoral joint of the knee caused by a degenerative, post-traumatic or inflammatory condition.
- Mobile version: osteoarthritis or avascular necrosis circumscribed to the internal tibiofemoral compartment of the knee.

5 - Implant selection

All of the implants that make up the X.NOV® knee prosthesis system must be used with implants from the same X.NOV® manufacturer.
The pre-operative planning is an essential step and is absolutely necessary to ensure the best clinical result. Comply with the surgical technique provided by X.NOV®.
X-ray photography allows the surgeon to plan the deformation angles to be transferred to the tibial instrumentation. Make sure to take the scale mentioned on the photographs into account.
It is absolutely necessary that you refer to the recommendations below in order to correctly choose implants that can be used together. An unsuitable choice of implants can lead to a reduction in the performance of the prosthesis, or even a reduction in its useful life.
It is essential to adhere strictly to the colour code on the device boxes. For mobile inserts to be combined with femoral components, devices with the same colour code (same colour) should always be selected. Any other assembly will have negative consequences on the performance of the device, including an increase in wear-and-tear and a decrease in the stability of the knee prosthesis.
KAPS® implants are available in versions with or without cement. It is strictly inadvisable to use the cementless versions:

- on fragile bones.
- on osteoporotic bones.
- in case of doubt regarding the primary stability of the components.
- if there is still a sensitive space between the osteoconductive coating and the bone.

The femoral component is available in 4 anatomical sizes (from F1 to F4) and in versions with or without cement.
The tibial base exists in 5 sizes (T1 to T5) and 4 versions: symmetrical mobile to be cemented, anatomical mobile to be cemented, anatomical fixed to be cemented and anatomical fixed without cement. The mobile version is only available in the cemented version.
The CHC 3.5 cancellous screw is available in three lengths (25, 35 and 45 mm) and is always to be combined with the anatomical fixed tibial base without cement.

Tibial base	Symmetrical		Anatomical	
	Cement	Cementless	Cement	Cementless
Mobile	✓		✓	
Fixed			✓	✓ + cancellous screw

The fixed polyethylene insert is available in 5 thicknesses (9, 10, 11, 12 and 14 mm) and in 5 sizes (T1 to T5) in combination with the tibial base.
The mobile polyethylene insert is available in 5 thicknesses (8, 9, 10, 11 and 12 mm) and in 4 sizes (F1 to F4) in combination with the femoral component.
The FULL-PE tibial plateau is available in 4 thicknesses (8, 10, 12 and 14 mm) and in 5 sizes (T1 to T5).

Thickness	PE TIBIAL INSERTS		
	FIXED	MOBILE	FULL-PE
H08		✓	✓
H09	✓	✓	
H10	✓	✓	✓
H11	✓	✓	
H12	✓	✓	✓
H14	✓		✓

Warning:

A size N mobile insert must only be mounted on a size N femoral component. A size N fixed insert must only be mounted on a size N tibial base.
It is recommended to place the mobile version of the KAPS® unicompartmental prosthesis in the internal compartment only.
Correction of a unicompartmental prosthesis with another unicompartmental prosthesis is discouraged, except for the isolated replacement of the tibial insert.

6 - Precautions for use

The user must have a thorough understanding of devices, instruments and surgical technique before the operation.
Inappropriate selection, positioning or securing of implants can lead to abnormal strain and reduce their useful life. Implantation must always be performed with the specific instruments provided by X.NOV® and in strict adherence with the surgical technique. The cementing technique is very important: comply with the surgical technique provided by X.NOV®. It is strictly forbidden to make any changes to the implant.
The operation must always be followed by appropriate care and rehabilitation.

Before the operation, check that:

- All the implants that make up the knee prosthesis are available and compatible with each other (check the characteristics indicated on the labelling and leaflet).
- The instruments to be used for the implantation are available and functional. Twisted or damaged instruments can cause improper positioning of the implants and defectiveness.
- Plan for implants that are at least one size larger and smaller than those determined before the procedure.
- It is recommended that the availability of a TKR be checked prior to the operation with a unicompartmental prosthesis to allow for a change in strategy during surgery if necessary.

General warnings:

- Only remove the implant from the package (in compliance with the usual sterile technique) once the correct size has been determined.
- Pay particular attention to all polished surfaces acting as friction surfaces, and to coated or treated surfaces.
- Never use implants that have been damaged, explanted, subjected to knocks or friction against a hard material.
- The dimensions of the blades used must comply with the X.NOV® specifications: thickness 1.27 mm, width 13 mm.
- Use of cemented versions is recommended in fragile and/or osteoporotic bones.
- Carefully check freedom of movement in order to detect incorrect fixations, instabilities or impingements; correct if necessary. Foreign particles can cause excessive wear-and-tear of the metal/plastic interface. The presence of osteophytes, cement residue and/or bone debris can cause the dislocation of the device or limit movements and make them painful.
- If cement is used, refer to the manufacturer's user's manual.

Warning:
Clean the surfaces well. They must be free from any foreign bodies prior to assembly of the various components and prior to the closure of the surgical site.

DO NOT IMPLANT:
- if the expiry date has been passed.
- if the packaging is damaged.

7 - Contraindications

- Arthritic damage to the patellofemoral compartment.
- Arthritic damage to the other tibiofemoral compartment.
- Excessive deformity of the patient's limb: genu varum or genu valgum > 10°.
- Mobile insert version: deformity of the limb in genu varum or genu valgum > 6°.
- Non-reducible flaccid or significant recurvatum.
- Flaccid greater than 15°.
- Disorder of the central pivot (ACL, PCL) and/or lateral ligaments (LLI, LLE).
- Severe muscle, nerve or vessel deficiencies affecting the limb in question.
- Bone fracture and/or major bone defects and/or poor bone quality, which may affect the stable implantation of the prosthesis.
- Any concomitant disease which might influence the implant function.
- Known allergy to implant components.
- Acute or chronic local or general infections (septicaemia, osteomyelitis...).
- Chronic joint instability (significant instability due to an advanced loss of osteo-cartilaginous structure or a lack of integrity of the lateral ligaments).
- Serious illnesses entailing a risk of dangerous postoperative complications.
- Rheumatoid polyarthritis and other forms of inflammatory joint disease.
- Inflammatory rheumatismal and microcrystalline joint damage.
- Severe osteoporosis.
- Osteomyelitis.
- Upper tibial osteotomy.
- Bone tumours.
- Poor skin coverage or persistent skin lesion on the knee joint.

Infection, septicaemia, osteomyelitis and a missing anterior cruciate ligament are absolute contraindications.

8 - Factors likely to compromise the success of the implantation

- Intense physical activities, risky activities or repetitive movements that expose the implant to excessive strain.
- Allergy to the materials of the implant.
- Medical history of infections and falls.
- BMI > 30 kg/m².
- Diabetes.
- Serious ligament instability and untreatable muscle contractures.
- Severe osteoporosis.
- History of osteotomy.
- Uncooperative or neurologically impaired patient.
- Mental or neuromuscular disorders that could constitute an unacceptable risk for the patient and that could cause post-operative complications.
- Insufficient bone quality.
- Direct contact of hydroxyapatite-coated implants with acidic local antibiotics.
- Metabolic diseases.
- Drug addiction and/or tendency to abuse drugs, medicinal products, alcohol or tobacco.
- Local bone tumours.
- Hypercorrection of the knee.
- Damage to implants during surgery.
- Interposition of a foreign body (blood, debris, shavings, etc.) between the implants.
- Pre-operative planning incorrect or not conducted.
- Incorrect selection of the size or type of implants used.
- Non-compliance with the surgical technique.
- Incorrect positioning of the implants.
- Use with an unsuitable instrument or device.
- Assembly of an insert with a femoral or tibial component of different sizes.

9 - Adverse effects

- Tibial collapse.
- Lengthening or shortening of the limb.
- Noises, squeaking.
- Anaphylactic shock.
- Loosening of one or more components.
- Dislocation due to trauma or wear of an UHMW PE tibial plateau.
- Pain.
- Soft tissue damage.
- Femoral or tibial supracondylar fracture.
- Stress fracture of the implant.
- Bone fracture during surgery.
- Infection.
- Instability of the implant.
- Infarction.
- Haematoma.
- Haemarthrosis.
- Hydrarthrosis.
- Radiological borders.
- Migration of the implant and loosening of the prosthesis.
- Metallosis.
- Osteolysis.
- Loss of mobility.
- Wound, poor healing.
- Revision.
- Reaction of tissue to debris from wear-and-tear.
- Reduced range of movement.
- Sensitivity reactions to the implanted materials.
- Rigidity of the knee.
- Sub-dislocation or dislocation of the mobile plateau.
- Cardiovascular and thrombo-embolic disorders: venous thrombosis, phlebitis and pulmonary or fatty embolism.
- Neuro-algo-dystrophic syndrome.
- Wear-and-tear of components.

And all of the traditional and known complications of prosthetic surgery.
Revision surgery may be necessary as a result of the inevitable wear and tear of the device, or any separation that may occur over time.

10 - Storage and handling

Implants must be stored unopened in their original packaging in a cool, dry place away from light. They must not be exposed to ionising radiation, extreme temperatures or particle contamination.
Implants are individually packed in double packaging. Make sure that the protective packaging is not damaged before removing the implant, as such damage can compromise the sterility of the implant.
Only remove the devices from the package immediately before use.
Disposal of the implant should be carried out in accordance with the method of waste disposal appropriate for the healthcare facility.
Warnings on the package as well as pictograms should be taken into account.

11 - Sterility

The implants are delivered sterile (gamma radiation at a minimum dose of 25 kGy).
They must be kept in their sealed original packaging until use.
Check the sterility expiry date on the label and check the protective packaging for signs of deterioration. If the protective packaging is damaged or the sterility expiry date exceeded, the implant should be sent back to X.NOV®. The rules for asepsis must be followed when removing the implant from the original packaging and for the entire duration of implantation.
Implants must not be re-sterilised or re-used, as this could lead to contamination and/or infection of the patient and damage to the device's properties.
Sterility is compromised if the packaging is damaged.

12 - Patient information

The surgeon must inform the patient of surgical risks, contraindications, possible adverse effects and potential interactions with other procedures or therapeutic or diagnostic devices (e.g.: MRI). He/she must talk to the patient about all the physical and psychological limitations inherent in the use of the implanted device. Particular attention must be paid to the patient's level of activity, weight, and need for regular medical check-ups.
The surgeon must warn the patient that the implant cannot be as robust, reliable and durable as normal and healthy bones. Therefore, the implant does not allow for the same levels of activity and constraints as those enabled by normal and healthy bones.

The surgeon must explain to the patient that the implant can break or be damaged following a fall or the practice of intense physical activity that could cause trauma or abnormal wear-and-tear to the implant.
If the patient has a job or activity that puts excessive strain on the implant (e.g. skiing, running, climbing, tennis, body building, etc.), the surgeon must warn the patient and inform him or her that the force exercised on the device and discontinuous and sudden movements that involve the prosthetic joint can cause the failure or reduction in the useful life of the implant or its revision.
Heavy temporary overload in case of falls, accidents or excessive overloads can cause the rupture of the implant, sometime even long after such an event.
The surgeon must tell the patient of the need to keep his/her weight under control and to have regular clinical and radiographic check-ups after the operation.

If the patient experiences pain, loss of mobility or other signs that would lead him or her to believe that the device is not or no longer performing as expected, he or she must consult the practitioner.

13. Compatibility and safety with MRI environment

Compatibility and safety with MRI environment of X.NOV® metallic implants have not been tested to evaluate risk of heating or migration.
Risks associated to MRI in presence of a metallic implant include migration and loosening, heating of the implant, heating and damage of the surrounding tissues, and artefactual distortion, specifically in the area surrounding the implant. Interpretation of MRI scans may be compromised or made impossible by such artefactual distortion.
Patients must be examined by a professional who is knowledgeable of the MRI equipment prior to undertaking MRI scanning.
Patients must inform the professional responsible for their examination of the presence of their prosthesis in case of MRI scanning.
Attention to any adverse effects that may be felt by the patient must be maintained throughout the MRI examination.

14. Instrumentation information

The instrumentation for insertion of the implants is supplied in a non-sterile state, unless the contrary is indicated on the packaging. The instrumentation is reusable and must be sterilised by the healthcare facility before use. See the leaflet and the instructions supplied with the instrumentation for information on decontamination and vapour sterilisation of the instrumentation.
Please contact X.NOV® for any additional information on the use of this device.

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